


TD\_Phosphate\_N1\_V2\_E

## Cellufine™ Phosphate

### Introduction

Cellufine™ Phosphate is the new product of Cellufine™ series. Cellufine™ Phosphate is a porosity bead although the Cellulose Phosphate known from ancient times is a fiber. Therefore, Cellufine™ Phosphate is the fast flow rate and high resolution compared with the conventional product.

Cellufine™ Phosphate and cellulose phosphate micrograph



Cellufine™ Phosphate



Conventional Cellulose Phosphate

### Comparison of resolution




Fig.1 Separation of mixed sample

Column Size: ID 1.1cm – Height 10cm

Flow rate: 2ml/min (126cm/h)

Buffer: 0.01M acetate buffer, pH4.8

Elution: 0 to 1 mol/L NaCl gradient

Sample

- ① ATP 0.7 mg/mL
- ② Ovalbumin 20 mg/mL
- ③ Lysozyme 8 mg/mL
- ④ Chytochrome C 10mg/mL

### Comparison of pore size





Fig.2 Cellufine Phosphate has a large pore compared with conventional Cellulose Phosphate.

### Fig.3 C.I.P. Stability




Cellufine Phosphate is stable to cleaning in place (C.I.P) by alkali.

## Cellufine Phosphate

Cellufine Phosphate is preferably applicable to purification of DNA binding protein.

This data was kindly provided by Dr. Svetlana Sedelnikova of University of Sheffield.



**Fig. Use of Cellufine Phosphate in Rus A D70N purification**

### Chromatography

Column: 1.6x10cm (20ml) packed with Cellufine Phosphate      Flow rate: 3ml/min( 90cm/h )

Sample: 7.5mg of Rus A D70N obtained after Heparin-Sepharose chromatography

Gradient: 200ml from 0.1 to 1.3M NaCl in 50mM tris-HCl pH 8.0

### SDS-PAGE

Gel: Novex 4-12%BT gel used with MES-SDS running buffer (Invitrogen)

Key for gel: 1 Cell free extract

2 Unbound material from Heparin-Sepharose column

3 RusA sample obtained from Heparin-Sepharose column

4-6 Fractions across the peak eluted from Cellufine Phosphate

7 RusA reference sample

8 Mark 12 MW standard (Invitrogen)

Fractions shown on lines 5 and six were combined to produce 3mg of the RusA D70N suitable for crystallization experiments (about 90% pure).

Ref.

Nucleic Acids Research, 2006, Vol. 00, No. 00 1–8

Rachel Macmaster, Svetlana Sedelnikova, Patrick J. Baker, Edward L. Bolt1, Robert G. Lloyd1 and John B. Rafferty

RusA Holliday junction resolves: DNA complex structure—insights into selectivity and specificity