

Super Edge Empty Mini Column

Versatility of empty mini column and performance evaluation of packed column- 2

JNC's empty mini columns can be packed with chromatography media from various manufacturers. In column performance evaluation-1, the quality of packed column was evaluated by theoretical plate number and symmetry factor. In this performance evaluation-2, we will introduce the data on the reliability of the adsorption capacity using the *Super Edge* empty mini column.

1 Compatibility; Comparison of adsorption capacities of commercially available packed columns and columns packed in *Super Edge* empty mini columns

Packed columns are sold by major chromatography resin manufacturers. The resins packed in these commercially available packed columns were purchased in bottles and packed in *Super Edge* empty mini columns, and the adsorption capacities of these were compared. If the comparison data show equivalence, the *Super Edge* empty mini column is compatible with the commercial packed column. It also shows that the adsorption performance of the resin of each company can be measured correctly.

Packing method

Packing was performed according to the manual using the packing tool enclosed in the starter kit. The manual can be obtained from *Super Edge* website.

<https://www.jnc-corp.co.jp/fine/se/index.html>

Material

Commercially available packed columns and bottled resins shown below were used.

<<Table1 AEX>>


| Manufacturer | Product abbreviation | Bead size μm | Base material |
|--------------|----------------------|-------------------------|---------------|
| T | DEAE(1) | 65 | Polymer |
| | Q(2) | 65 | Polymer |
| | Q(3) | 75 | Polymer |
| | QAE(4) | 100 | Polymer |
| G | DEAE(5) | 90 | Agarose |
| | Q(6) | 90 | Agarose |
| B | Q(7) | 50 | Polymer |

<<Table2 CEX>>

| Manufacturer | Product abbreviation | Bead size μm | Base material |
|--------------|----------------------|-------------------------|---------------|
| T | CM(8) | 75 | Polymer |
| | S(9) | 75 | Polymer |
| | SP(10) | 75 | Polymer |
| | SP(11) | 75 | Polymer |
| G | S(12) | 90 | Agarose |

Chromatography

The dynamic adsorption capacity was measured by performing chromatography under the conditions of Table 3 and 2. A typical chromatogram is shown below. Loading was performed until 10% of the protein concentration was detected, followed by washing and then elution to determine the amount of adsorption and recovery.



<< Table 3 Chromatography condition >>

| Step | CV(column volume) ml |
|-----------------------|-------------------------|
| 1. Equilibration | 5 |
| 2. Sample application | Up to 10% break through |
| 3. Column wash | 10 |
| 4. Elution1 (0→100%) | 20 |
| 5. Elution 2(100%) | 10 |
| 6. Column wash | 10 |
| 7. CIP | 10 0.1M NaOH |
| 8. Equilibration | 20 |

※Flow velocity: Column residence time 1
1ml column 1ml/ml, 5ml column 5ml/ml

<< Table 4 Protein and buffer condition >>

| Mode | Protein | Buffer A | Buffer-B |
|------|-----------------|---------------------------------|--------------------|
| AEX | BSA 5mg/ml | 50mM Tris-HCl, pH8.5 | Buffer-A + 1M NaCl |
| CEX | Lysozyme 5mg/ml | 10mM Na-sodium phosphate, pH7.0 | Buffer-A + 1M NaCl |

JNC CORPORATION

Result

Table 5 Comparison of 10% DBC (mg / ml) and recovery rate of *Super Edge* mini column and commercial packed column

| | 10% DBC | | Recovery (%) | |
|---------|----------|------------|--------------|------------|
| | JNC | Commercial | JNC | Commercial |
| AEX | BSA | | | |
| DEAE(1) | 25.2 | 21.4 | 105 | 111 |
| Q(2) | 112.7 | 74.5 | 98 | 98 |
| Q(3) | 124.8 | 134.5 | 100 | 98 |
| QAE(4) | 26.9 | 19.5 | 103 | 118 |
| DEAE(5) | 166.2 | 159.6 | 101 | 101 |
| Q(6) | 120.2 | 111.9 | 97 | 95 |
| CEX | Lysozyme | | | |
| CM(8) | 151.9 | 132.8 | 101 | 101 |
| S(9) | 155.0 | 142.9 | 100 | 100 |
| SP(10) | 78.4 | 65.6 | 110 | 109 |

※JNC : *Super Edge* mini column,
 ※Commercial: Packed column purchased from a resin manufacturer.




Fig1. Regression analysis of 10% DBC (mg / ml) of *Super Edge* mini column and commercial column

Comparative data are shown in Table 5 and Fig 1. Q (2) and Q (3) seem to be different in DBC, but the result of regression analysis showed that the intercept was about 1 and the correlation coefficient was close to 1, indicating compatibility.

JNC CORPORATION

2 1 ml and 5 ml scalability of *Super Edge* empty mini column

There are two types of JNC's *Super Edge* empty mini columns, 1 ml and 5 ml, and it was investigated whether the adsorption capacity could be measured equally on these scales.

<< Table 6 Column dimension >>

| Column volume ml | I.D. mm | Length mm |
|------------------|---------|-----------|
| 1 | 6.7 | 30 |
| 5 | 14.6 | |

Material

The resins in Tables 1 and 2 and the following ion exchangers were tested.

<< Table 6 Resin added to comparative experiment >>

| Manufacturer | Product abbreviation | Bead size μm | Base material |
|-------------------|----------------------|-------------------------|---------------|
| JNC Cellufine® | A-500 | 90 | Cellulose |
| | MAX DEAE | 90 | Cellulose |
| | MAX Q-h | 90 | Cellulose |
| | MAX Q-r | 90 | Cellulose |
| | MAX S-h | 90 | Cellulose |
| | MAX S-r | 90 | Cellulose |
| T | S(13) | 100 | Polymer |
| G | S(14) | 90 | Agarose |

Chromatography

The procedure was similar to that described above.

Result

<< Table7 Comparison of DBC with 1ml and 5ml *Super Edge* mini column>>

| | | 10%DBC(mg/ml) | |
|-----------|----------|---------------|------------|
| | | 1ml Column | 5ml Column |
| Cellufine | A-500 | 18 | 18 |
| | MAX DEAE | 154 | 166 |
| | MAX Q-h | 174 | 202 |
| | MAX Q-r | 131 | 137 |
| | MAX S-h | 165 | 192 |
| | MAX S-r | 126 | 122 |
| | T | DEAE(1) | 25 |
| Q(2) | | 113 | 129 |
| Q(3) | | 125 | 119 |
| QAE(4) | | 27 | 27 |
| CM(8) | | 152 | 140 |
| S(9) | | 155 | 142 |
| SP(10) | | 78 | 78 |
| S(13) | | 30 | 29 |
| G | DEAE(5) | 120 | 119 |
| | Q(6) | 166 | 164 |
| | S(14) | 151 | 139 |
| B | Q(7) | 36 | 33 |




Fig2 Regression analysis of DBC on 1ml and 5ml *Super Edge* mini columns

As shown in Table 7 and Figure 2, the *Super Edge* empty mini columns had similar adsorption capacities measured at 1 ml and 5 ml.

JNC CORPORATION

3 CIP cycle and stability

We investigated the relationship between the number of times a packed column made with *Super Edge* mini-column was used and its performance. The column was loaded with protein → recovery → CIP → equilibration was repeated, and changes in theoretical plate number, symmetry factor, adsorption amount, and recovery rate were investigated.

Chromatography

Under the conditions of CEX in Table 4, the cycles 1 to 8 of Table 3 were repeated, and the theoretical plate number and symmetry factor *As* were examined for each cycle.

Column

Cellufine® MAX S-h was packed in 1 ml and 5 ml respectively.

Measurement

Measurements were performed once for each cycle on 3 columns of 1 ml and 5 ml, and the average was calculated from the values of 3 columns.

Result

<< Table 8. 1ML COLUMN / CIP cycle and column performance >>

| Cycle | N/m | | As | | 10%DBC(mg/ml) | | Recovery | |
|-------|---------|------|---------|------|---------------|------|----------|------|
| | Average | S.D. | Average | S.D. | Average | S.D. | Average | S.D. |
| 0 | 5,738 | 421 | 0.93 | 0.07 | | | | |
| 1 | 5,488 | 246 | 0.87 | 0.01 | 172.6 | 4.0 | 98 | 0.5 |
| 2 | 5,366 | 86 | 0.86 | 0.06 | 176.2 | 4.0 | 96 | 0.8 |
| 3 | 5,593 | 210 | 0.85 | 0.05 | 165.6 | 3.5 | 103 | 1.6 |
| 4 | 5,340 | 180 | 0.86 | 0.02 | 171.3 | 2.1 | 100 | 0.2 |
| 5 | 5,542 | 226 | 0.86 | 0.02 | 168.3 | 3.7 | 102 | 0.5 |
| 6 | 5,607 | 169 | 0.86 | 0.01 | 172.0 | 4.1 | 98 | 0.7 |
| 7 | 5,523 | 100 | 0.90 | 0.01 | 171.4 | 4.0 | 99 | 0.4 |
| 8 | 5,636 | 130 | 0.87 | 0.02 | 166.5 | 4.8 | 105 | 0.4 |
| 9 | 5,490 | 30 | 0.88 | 0.03 | 167.4 | 3.2 | 102 | 0.7 |
| 10 | 5,639 | 118 | 0.91 | 0.05 | 174.7 | 3.6 | 96 | 0.4 |

<< Table 9. 5ML COLUMN / CIP cycle and column performance >>

| Cycle | N/m | | As | | 10%DBC(mg/ml) | | Recovery | |
|-------|---------|------|---------|------|---------------|------|----------|------|
| | Average | S.D. | Average | S.D. | Average | S.D. | Average | S.D. |
| 0 | 6,849 | 176 | 1.22 | 0.01 | | | | |
| 1 | 6,741 | 279 | 1.23 | 0.04 | 174.3 | 1.6 | 99 | 1.3 |
| 2 | 6,956 | 327 | 1.21 | 0.02 | 172.5 | 2.2 | 103 | 2.2 |
| 3 | 7,049 | 352 | 1.25 | 0.02 | 173.0 | 1.6 | 102 | 3.1 |
| 4 | 6,634 | 547 | 1.29 | 0.04 | 171.4 | 1.7 | 103 | 1.3 |
| 5 | 7,223 | 294 | 1.23 | 0.02 | 171.8 | 1.8 | 102 | 1.8 |
| 6 | 7,252 | 230 | 1.24 | 0.02 | 175.5 | 1.7 | 101 | 1.8 |
| 7 | 7,200 | 217 | 1.24 | 0.03 | 173.5 | 1.7 | 100 | 0.6 |
| 8 | 7,264 | 305 | 1.25 | 0.04 | 173.3 | 1.7 | 101 | 1.2 |
| 9 | 7,281 | 238 | 1.24 | 0.01 | 168.5 | 1.6 | 104 | 0.2 |
| 10 | 7,297 | 198 | 1.25 | 0.03 | 176.4 | 1.7 | 100 | 1.4 |




Fig3 CIP cycle and column performance

Tables 8 and 9 show the average value and standard deviation of the results measured for each cycle for 3 columns of 1 ml and 5 ml, respectively. Fig3 is a graph of Tables 8 and 9.

There was no change in the theoretical plate number, symmetry factor, and DBC and recovery rate for each CIP cycle, and no deterioration was observed up to 10 cycles. As far as the tendency of the theoretical plate number and symmetry factor in Fig. 3 is concerned, it seems that there is no sudden deterioration even after 10 cycles. However, care must be taken as DBC is affected by the durability of the packed resin. Please refer to the instruction manual of the resin you actually use and set the CIP conditions before use.

Conclusion

- The commercially available packed column and *Super Edge* mini column have the same recovery rate as DBC, and it is possible to purchase a commercially available resin bottle and packed it into the *Super Edge* mini column for performance evaluation.
- Since the DBC and recovery rates of 1 ml and 5 ml of *Super Edge* mini-column are the same, it was confirmed that scale-up can be performed easily.
- The packed column of the *Super Edge* mini column does not show any deterioration in performance after 10 repeated CIP cycles, and can perform several times chromatography experiments.