


What is Cellufine™ MLP?


Cellulose resin MLP is a highly crosslinked cellulose particle that has exceeding 1 μm continuous pore structure.

Therefore, the intraparticle surface of the resin can be used as an adsorption site for biomacromolecules such as large virus and virus-like particles.



- ✓ Improved mass transfer
- ✓ High binding capacity for large biomolecules



Green fluorescent nanoparticles (100 nm) can diffuse into the intraparticle area of Cellufine™ MLP.



Target molecule size



We welcome requests for free samples and collaborative research!

10 mL bottle



and/or

Prepacked column
(1 mL, 5 mL)



Let's create the future
together. Reach us via
the QR code on back.






Virus purification

High binding capacity of >100 nm bioparticles


Cellufine™ MLP DexS

Binding capacity of influenza virus



Column vol : ID 5 mm × 1.5 cm L (CV = 0.29 mL)
Flow rate : 0.5 mL/min
Load sample : Influenza A virus (H1N1) from MDCK cells

SARS-CoV-2 Purification



Column vol : ID 5 mm × 1.5 cm L (CV = 0.29 mL)
Flow rate : 0.5 mL/min
Load sample : SARS-CoV-2 XBB




AAV purification

High performance for AAV full/empty capsids separation

Cellufine™ MLP AEX

| | Full particle ratio* |
|------------------------|----------------------|
| Load (AFF purified) | 15.9% |
| Monolithic Column | 42.8% |
| MLP AEX for AAV | 70.6% |

* Full particle ratio was determined using a mass photometry



Column vol : 1.0 mL

Flow rate : 1.0 mL/min (RT 1min)

Load sample : AFF purified AAV2: 1.06×10^{12} vg

Equilibration buffer : 50 mM Tris-HCl, pH 9.0, 2 mM MgCl₂

Elution buffer (B) : 50 mM Tris-HCl, 150 mM NaCl, 2 mM MgCl₂, pH 9.0



mAb purification

Efficient host cell proteins removal under high loading condition

Cellufine™ MLP MMC

| | mAb purity in feed (%) | HCPs in feed (ppm) | mAb purity in FT pool (%) | mAb yield (%) | HCPs in FT pool (ppm) |
|----------------|------------------------|--------------------|---------------------------|---------------|-----------------------|
| Capto adhere | 97.2 | 1718 | 97.6 | 95.6 | 513 |
| MAX IB | 96.7 | 2229 | 97.8 | 96.3 | 59 |
| MLP MMC | 97.2 | 1718 | 97.7 | 99.4 | 8 |

Column vol : 0.3 mL

Flow rate : 0.075 mL/min

Load sample : Partially purified mAb solution (pH 7.0, 6.0 mS/cm), 1010 mg mAb/1 mL resin