

JNC CORPORATION

2-1, Otemachi 2-chome, Chiyoda-ku, Tokyo, 100-8105, Japan

Phone: +81-3-3243-6150 / Facsimile: +81-3-3243-6219

March 6, 2018

Dear Cellufine Customer :

This letter is to inform you of a minor change to 5 mL of Cellufine mini-column products that you have purchased. If you have questions regarding this change, please contact the JNC staff referenced below.

Detailed Information of the PCN (Product Change Notification)

PCN Title: 5 mL of Cellufine mini-column
Description of the Change: Design, connector parts and size of appearance


<Change of Product>

Current product



13 mm I.D. x 44 mm L (5 mL)

New product



15 mm I.D. x 30 mm L (5 mL)

Impact and Recommended Action:

Any Cellufine resins do not change.

Effective Date:

April 1, 2018

Who to Contact with Questions:

If you have questions and concerns, please contact at

cellufine@jnc-corp.co.jp or please visit our web site

<http://www.jnc-corp.co.jp/fine/en/cellufine/index.html>

Best regards,

Toshiharu Iwata



Technical Manager, Life Chemicals Launch Office