

JNC CORPORATION

2-1, Otemachi 2-chome, Chiyoda-ku, Tokyo, 100-8105, Japan

Phone: +81-3-3243-6150 / Facsimile: +81-3-3243-6219

February 15, 2017

Dear Cellufine Customer:

This letter is to inform you of a minor change to 1mL of Cellufine mini-column products that you have purchased. If you have questions regarding this change, please contact the JNC staff referenced below.

Detailed Information of the PCN (Product Change Notification)

PCN Title: 1 mL of Cellufine mini-column
Description of the Change: Design, connector parts and size of appearance

<Change of Product>

Current product



9.0 mm x 18 mm (1 mL)

New product



6.7 mm x 30 mm (1 mL)

Impact and Recommended Action: Any Cellufine resins do not change.
Effective Date: March 1, 2017
Who to Contact with Questions: If you have questions and concerns, please contact at cellufine@jnc-corp.co.jp or please visit our web site <http://www.jnc-corp.co.jp/fine/en/cellufine/index.html>

Best regards,

A handwritten signature in black ink, appearing to read 'T. Iwata'.

Toshiharu Iwata
Technical Manager, Life Chemicals Launch Office