

JNC CORPORATION

Shin-Otemachi Bldg.,
2-1, Otemachi 2-Chome
Chiyoda-ku, Tokyo 100-8105 JAPAN
Life Chemicals Division

Phone : +81-3-3243-6225 Facsimile : +81-3-3243-6219 E-mail: mcmp@jnc-corp.co.jp

February 1, 2026

Dear Valued Customer,

Notification of Shelf-Life Extension and Change in Storage Temperature for MC-Media Pad SA

This letter is to notify you of the shelf-life extension and the change in storage temperature for the following product:

Product name:

MC-Media Pad SA

Details of the Change:

- **Reason for Change:** The shelf-life extension is based on successful stability testing, which has confirmed that product quality is maintained throughout the new shelf-life period. In addition, the storage temperature will be changed to further ensure the continued quality of the product.
- **Package Change:** Due to the change in storage temperature, the package will also be changed.
- **Effective Date:** This change will apply to production lots manufactured on or after March 1, 2026.

	Before	After
Shelf-Life	1 year after production	3 years after production
Storage Temp.	2-15°C	2-8°C
Aluminum Bag	<p>MC-Media Pad™ サニ太くん®/Sanita-kun</p> <p>SA 03 Staphylococcus aureus 黄色ブドウ球菌用 金黄色葡萄球菌用 黄緑芽生球菌用</p> <p>Incubation condition 35°C 24h</p> <p>2-15°C 25pcs</p>	<p>MC-Media Pad™ サニ太くん®/Sanita-kun</p> <p>SA 03 Staphylococcus aureus 黄色ブドウ球菌用 金黄色葡萄球菌用 黄緑芽生球菌用</p> <p>Incubation condition 35°C 24h</p> <p>2-8°C 25pcs</p>
Paper Box	<p>MC-Media Pad™ サニ太くん®/Sanita-kun</p> <p>SA 03 Staphylococcus aureus 黄色ブドウ球菌用 金黄色葡萄球菌用 黄緑芽生球菌用</p> <p>Incubation conditions 35°C 24h</p> <p>2-15°C 25pcs x4</p>	<p>MC-Media Pad™ サニ太くん®/Sanita-kun</p> <p>SA 03 Staphylococcus aureus 黄色ブドウ球菌用 金黄色葡萄球菌用 黄緑芽生球菌用</p> <p>Incubation conditions 35°C 24h</p> <p>2-8°C 25pcs x4</p>
Carton (All Products)	<p>JNC CORPORATION MADE in JAPAN</p> <p>2-15°C 1000pcs</p> <p>MC-Media Pad™ サニ太くん®/Sanita-kun</p>	<p>JNC CORPORATION MADE in JAPAN</p> <p>冷蔵 1000pcs</p> <p>MC-Media Pad™ サニ太くん®/Sanita-kun</p>

JNC CORPORATION

Shin-Otemachi Bldg.,
2-1, Otemachi 2-Chome

Chiyoda-ku, Tokyo 100-8105 JAPAN

Life Chemicals Division

Phone : +81-3-3243-6225 Facsimile : +81-3-3243-6219 E-mail: mcmp@jnc-corp.co.jp

Please Note:

- To further ensure product quality, the storage temperature range will be changed to 2–8°C, and packaging will be updated accordingly. Eventually, all products will use 2–8°C storage; for now, only cartons will be changed.
- Although previously shipped MC-Media Pad SA products may have the old conditions printed on the packaging, products that have been stored at 2–8°C may be used according to the new conditions described in this notification (shelf-life: 3 years after production / storage temperature: 2–8°C).
- Please note that products stored under the old conditions (storage temperature: 2–15°C) should remain subject to the previous shelf-life -1 year after production.
- This change does not affect product quality or performance in any way.

If you have any questions or require further information, please do not hesitate to contact us.
Thank you for your continued support of the MC-Media Pad series.

Sincerely,
Life Chemicals Division
JNC CORPORATION